




Fitzgerald

Research Products Catalog

2007 - 2008



Retroviral Products

Recombinant Proteins

Monoclonal Antibodies

Polyclonal Antibodies

HIV-1 Recombinant Proteins

Antigen	Catalog #	Expression System	Purity as estimated by SDS-PAGE Gel	Label	Qty	Price
gp120 IIIB	30-12001-2	Baculovirus	>90% Purity	Unlabeled	100µg	\$504
gp120 IIIB	65-012002	Baculovirus	>90% Purity	HRP labeled	50µg	\$610
gp120 IIIB	65-012003	Baculovirus	>90% Purity	FITC labeled	50µg	\$610
gp41 IIIB (amino acids 1-241)	30-036001	E.coli	>98% Purity	Unlabeled	100µg	\$461
gp41 IIIB (amino acids 51-131)	30-039001	E.coli	>90% Purity	Unlabeled	100µg	\$495
gp41 MN (amino acids 52-258)	30-040001	E.coli	>95% Purity	Unlabeled	100µg	\$461
gp41 MN (amino acids 52-258)	65-040002	E.coli	>95% Purity	HRP labeled	50µg	\$524
gp41 MN (amino acids 52-258)	65-040003	E.coli	>95% Purity	FITC labeled	50µg	\$524
gp41 MN (amino acids 52-258)	65-040005	E.coli	>95% Purity	Biotin labeled	50µg	\$524
gp160 IIIB soluble	65-024102	Baculovirus	>90% Purity	HRP labeled	50µg	\$525
p24 IIIB (amino acids 1-258)	65-014002	Baculovirus	>90% Purity	HRP labeled	50µg	\$525
p24 IIIB (amino acids 1-258)	65-014003	Baculovirus	>90% Purity	FITC labeled	50µg	\$525
p24 IIIB (amino acids 1-258)	65-014005	Baculovirus	>90% Purity	Biotin labeled	50µg	\$525
p24 HxB2 (amino acids 1-189)	30-014101	E.coli	>90% Purity	Unlabeled	100µg	\$461
p24 HxB2 (amino acids 1-189)	65-014102	E.coli	>90% Purity	HRP labeled	50µg	\$525
p24 HxB2 (amino acids 1-189)	65-014105	E.coli	>90% Purity	Biotin labeled	50µg	\$525
p17 IIIB (amino acids 2-120)	65-048003	E.coli	>90% Purity	FITC labeled	50µg	\$576
p17 IIIB (amino acids 4-120)	30-049001	E.coli	>90% Purity	Unlabeled	100µg	\$461
Protease IIIIB	30-031001	E.coli	>90% Purity	Unlabeled	100µg	\$461
Protease IIIIB	65-031002	E.coli	>90% Purity	HRP labeled	50µg	\$525
Protease IIIIB	65-031005	E.coli	>90% Purity	Biotin labeled	50µg	\$525
RT IIIIB	30-026001	E.coli	>90% Purity	Unlabeled	20µg	\$694
tat HxB2	30-017001	E.coli	>90% Purity	Unlabeled	100µg	\$1055
tat HxB2	65-017002	E.coli	>90% Purity	HRP labeled	50µg	\$1161
tat HxB2	65-017003	E.coli	>90% Purity	FITC labeled	50µg	\$1161
tat HxB2	65-017005	E.coli	>90% Purity	Biotin labeled	50µg	\$1161
rev HxB2/HxB3	30-016001	E.coli	>90% Purity	Unlabeled	100µg	\$694
rev HxB2/HxB3	65-016003	E.coli	>90% Purity	FITC labeled	50µg	\$800
rev HxB2/HxB3	65-016005	E.coli	>90% Purity	Biotin labeled	50µg	\$800
nef ELI	30-025001	E.coli	>90% Purity	Unlabeled	100µg	\$694
nef ELI	65-025003	E.coli	>90% Purity	FITC labeled	50µg	\$737
nef ELI	65-025005	E.coli	>90% Purity	Biotin labeled	50µg	\$737
VIF IIIIB	30-028001	E.coli	>90% Purity	Unlabeled	100µg	\$595

HIV-2 Recombinant Proteins

Antigen	Catalog #	Expression System	Purity as estimated by SDS-PAGE Gel	Label	Qty	Price
gp105 ROD	30-015001	Baculovirus	>90% Purity	Unlabeled	100µg	\$843
gp105 ROD	65-015003	Baculovirus	>90% Purity	FITC labeled	50µg	\$949
gp105 ROD	65-015005	Baculovirus	>90% Purity	Biotin labeled	50µg	\$949
p26 ROD	30-018001	E.coli	>90% Purity	Unlabeled	100µg	\$694
p26 ROD	65-018002	E.coli	>90% Purity	HRP labeled	50µg	\$737
p26 ROD	65-018003	E.coli	>90% Purity	FITC labeled	50µg	\$737
gp36 ROD (amino acids 43-148)	30-034001	E.coli	>90% Purity	Unlabeled	100µg	\$410
gp36 ROD (amino acids 43-148)	65-034002	E.coli	>90% Purity	HRP labeled	50µg	\$495
gp36 ROD (amino acids 43-148)	65-034003	E.coli	>90% Purity	FITC labeled	50µg	\$495
rev ROD	30-019001	E.coli	>90% Purity	Unlabeled	100µg	\$795
RT ROD (amino acids 43-148)	30-027001	E.coli	>90% Purity	Unlabeled	20µg	\$694

SIV Recombinant Proteins

Antigen	Catalog #	Expression System	Purity as estimated by SDS-PAGE Gel	Label	Qty	Price
gp120 MAC251	30-033001	Baculovirus	>90% Purity	Unlabeled	100µg	\$912
gp120 MAC251	65-033005	Baculovirus	>90% Purity	Biotin labeled	50µg	\$1018
p28 MAC251 (amino acids 138-387)	30-022001	E.coli	>90% Purity	Unlabeled	100µg	\$694
p28 MAC251 (amino acids 138-387)	65-022002	E.coli	>90% Purity	HRP labeled	50µg	\$737
p28 MAC251 (amino acids 138-387)	65-022005	E.coli	>90% Purity	Biotin labeled	50µg	\$737
rev MAC251	30-021001	E.coli	>90% Purity	Unlabeled	100µg	\$694

HTLV-I/HTLV-II Recombinant Proteins

Antigen	Catalog #	Expression System	Purity as estimated by SDS-PAGE Gel	Label	Qty	Price
gp46 (amino acids 16-312)	30-041001	E.coli	>90% Purity	Unlabeled	100µg	\$560
gp46 (amino acids 26-215)	30-042001	E.coli	>90% Purity	Unlabeled	100µg	\$567
gp46 (amino acids 217-312)	30-043001	E.coli	>90% Purity	Unlabeled	100µg	\$567
gp21 (amino acids 14-121)	30-044001	E.coli	>90% Purity	Unlabeled	100µg	\$567
gp21 (amino acids 14-121)	65-044002	E.coli	>90% Purity	HRP labeled	50µg	\$631
gp21 (amino acids 14-121)	65-044003	E.coli	>90% Purity	FITC labeled	50µg	\$631
gp21 HTLV-II	30-046001	E.coli	>90% Purity	Unlabeled	100µg	\$630

Hepatitis B Virus Recombinant Proteins

Antigen	Catalog #	Expression System	Purity as estimated by SDS-PAGE Gel	Label	Qty	Price
Hepatitis B Core Protein (amino acids 1-144)	30-050001	E.coli	>95% Purity	Unlabeled	100µg	\$567

HIV-1 IIIB Monoclonal Antibodies

Antibody	Catalog #	Isotype	Epitope mapped to:	Label	Qty	Price
Mouse anti-gp120 Clone: α 70	10-001001	IgG ₁	gp120 V α loop	Unlabeled	100 μ g	\$223
Mouse anti-gp120 Clone: α 70	61-001002	IgG ₁	gp120 V α loop	HRP labeled	50 μ g	\$329
Mouse anti-gp120 Clone: α 70	61-001003	IgG ₁	gp120 V α loop	FITC labeled	50 μ g	\$350
Mouse anti-gp120 Clone: α 70	61-001005	IgG ₁	gp120 V α loop	Biotin labeled	50 μ g	\$345
Mouse anti-gp41 Clone: 1577	10-001201	IgG ₁	gp41 (733-750)	Unlabeled	100 μ g	\$235
Mouse anti-gp41 Clone: 1577	61-001202	IgG ₁	gp41 (733-750)	HRP labeled	50 μ g	\$313
Mouse anti-gp41 Clone: 1577	61-001203	IgG ₁	gp41 (733-750)	FITC labeled	50 μ g	\$313
Mouse anti-gp41 Clone: 1577	61-001205	IgG ₁		Biotin labeled	50 μ g	\$222
Mouse anti-gp41 Clone: #D1	61-002302	IgG ₁		HRP labeled	50 μ g	\$222
Mouse anti-gp41 Clone: #D1	61-002303	IgG ₁		FITC labeled	50 μ g	\$195
Mouse anti-gp41 Clone: #D1	61-002305	IgG ₁		Biotin labeled	50 μ g	\$285
Mouse anti-p24 Clone:26C10C10	61-004002	IgG _{2b}	p25 (154-172)	HRP labeled	50 μ g	\$313
Mouse anti-p24 Clone:26C10C10	61-004003	IgG _{2b}	p25 (154-172)	FITC labeled	50 μ g	\$313
Mouse anti-p24 Clone:26C10C10	61-004005	IgG _{2b}	p25 (154-172)	Biotin labeled	50 μ g	\$313
Mouse anti-P24 Clone: #18	10-004101	IgG	p25 (148-162)	Unlabeled	100 μ g	\$225
Mouse anti-P24 Clone: #18	61-004102	IgG	p25 (148-162)	HRP labeled	50 μ g	\$295
Mouse anti-P24 Clone: #18	61-004103	IgG	p25 (148-162)	FITC labeled	50 μ g	\$295
Mouse anti-P24 Clone: #18	61-004105	IgG	p25 (148-162)	Biotin labeled	50 μ g	\$295
Mouse anti-rev Clone: RMAb6	10-006001	IgG ₁	rev (96-110)	Unlabeled	100 μ g	\$345
Mouse anti-rev Clone: RMAb6	61-006002	IgG ₁	rev (96-110)	HRP labeled	50 μ g	\$419
Mouse anti-rev Clone: RMAb6	61-006005	IgG ₁	rev (96-110)	Biotin labeled	50 μ g	\$419
Mouse anti-rev Clone: 9G2G4D6E8	10-006101	IgG _{2b}	rev (70-84)	Unlabeled	100 μ g	\$345
Mouse anti-rev Clone: 9G2G4D6E8	61-006103	IgG _{2b}	rev (70-84)	FITC labeled	50 μ g	\$419
Mouse anti-rev Clone: 9G2G4D6E8	61-006105	IgG _{2b}	rev (70-84)	Biotin labeled	50 μ g	\$419
Mouse anti-tat Clone: tMAB-B	10-007001	IgG ₁	tat (6-14)	Unlabeled	100 μ g	\$345
Mouse anti-tat Clone: tMAB-B	61-007002	IgG ₁	tat (6-14)	HRP labeled	50 μ g	\$419
Mouse anti-tat Clone: tMAB-B	61-007003	IgG ₁	tat (6-14)	FITC labeled	50 μ g	\$419
Mouse anti-tat Clone: tMAB-B	61-007005	IgG ₁	tat (6-14)	Biotin labeled	50 μ g	\$419
Mouse anti-nef N-terminus Clone: 2H12	10-002001	IgG ₁	nef (54-67)	Unlabeled	100 μ g	\$223
Mouse anti-nef N-terminus Clone: 2H12	61-002003	IgG ₁	nef (54-67)	FITC labeled	50 μ g	\$313
Mouse anti-nef N-terminus Clone: 2H12	61-002005	IgG ₁	nef (54-67)	Biotin labeled	50 μ g	\$313
Mouse anti-nef C-terminus Clone: 3F2	10-002101	IgG ₁		Unlabeled	100 μ g	\$313
Mouse anti-nef C-terminus Clone: 3F2	61-002103	IgG ₁		FITC labeled	50 μ g	\$313
Mouse anti-nef C-terminus Clone: 3F2	61-002105	IgG ₁		Biotin labeled	50 μ g	\$223
Mouse anti-RT Clone: 8C4D7	10-009001	IgG ₁		Unlabeled	100 μ g	\$223
Mouse anti-RT Clone: 8C4D7	61-009002	IgG ₁		HRP labeled	50 μ g	\$313
Mouse anti-RT Clone: 8C4D7	61-009003	IgG ₁		FITC labeled	50 μ g	\$313
Mouse anti-RT Clone: 8C4D7	61-009005	IgG ₁		Biotin labeled	50 μ g	\$313
Mouse anti-Integrase Clone: #1	10-001301	IgG ₁		Unlabeled	100 μ g	\$223
Mouse anti-Integrase Clone: #1	61-001303	IgG ₁		FITC labeled	50 μ g	\$313

HIV-2 ROD Monoclonal Antibodies

Antibody	Catalog #	Isotype	Recognises:	Label	Qty	Price
Mouse anti-gp 105 Clone: #21	10-008001	IgG ₁	gp105 (Catalog #015001)	Unlabeled	100µg	\$223
Mouse anti-gp 105 Clone: #21	61-008002	IgG ₁	gp105 (Catalog #015001)	HRP labeled	50µg	\$313
Mouse anti-gp 105 Clone: #21	61-008003	IgG ₁	gp105 (Catalog #015001)	FITC labeled	50µg	\$313
Mouse anti-gp 105 Clone: #21	61-008005	IgG ₁	gp105 (Catalog #015001)	Biotin labeled	50µg	\$313

SIV MAC251 Monoclonal Antibodies

Antibody	Catalog #	Isotype	Recognises:	Label	Qty	Price
Mouse anti-p28 Clone: 3F7	10-002201	IgG ₁	SIV _{mac251} p28 (Catalog #022001)	Unlabeled	100µg	\$215
Mouse anti-p28 Clone: 3F7	10-02201-1	IgG ₁	SIV _{mac251} p28 (Catalog #022001)	Unlabeled	100µg	\$450
Mouse anti-p28 Clone: 3F7	10-02201-2	IgG ₁	SIV _{mac251} p28 (Catalog #022001)	Unlabeled	100µg	\$350
Mouse anti-p28 Clone: 3F7	61-002202	IgG ₁	SIV _{mac251} p28 (Catalog #022001)	HRP labeled	50µg	\$295
Mouse anti-p28 Clone: 3F7	61-002203	IgG ₁	SIV _{mac251} p28 (Catalog #022001)	FITC labeled	50µg	\$295
Mouse anti-p28 Clone: 3F7	61-002205	IgG ₁	SIV _{mac251} p28 (Catalog #022001)	Biotin labeled	50µg	\$295

HTLV-I Monoclonal Antibodies

Antibody	Catalog #	Isotype	Reacts with:	Label	Qty	Price
Mouse anti-tax	10-810001	IgG ₁	Recombinant tax (HTLV) in West. Blot	Unlabeled	100µg	\$295
Mouse anti-tax	61-810003	IgG ₁	Recombinant tax (HTLV) in West. Blot	FITC labeled	50µg	\$370

T-Cell Receptor Monoclonal Antibodies

Antibody	Catalog #	Isotype	Remarks:	Label	Qty	Price
Mouse anti-CD4 Clone: #45	10-003001	IgG ₁	Recognises CD4 (Catalog #013001)	Unlabeled	100µg	\$210
Mouse anti-CD4 Clone: #45	61-003002	IgG ₁	Recognises CD4 (Catalog #013001)	HRP labeled	50µg	\$310
Mouse anti-CD4 Clone: #45	61-003003	IgG ₁	Recognises CD4 (Catalog #013001)	FITC labeled	50µg	\$310
Mouse anti-CD4 Clone: #45	61-003005	IgG ₁	Recognises CD4 (Catalog #013001)	Biotin labeled	50µg	\$312
Mouse anti-CD4 Clone: #4	10-003101	IgG ₁	Leu-3A equivalent	Unlabeled	100µg	\$210
Mouse anti-CD4 Clone: #4	61-003102	IgG ₁	Leu-3A equivalent	HRP labeled	50µg	\$310
Mouse anti-CD4 Clone: #4	61-003103	IgG ₁	Leu-3A equivalent	FITC labeled	50µg	\$310
Mouse anti-CD4 Clone: #4	61-003105	IgG ₁	Leu-3A equivalent	Biotin labeled	50µg	\$310
Mouse anti T-cell receptor CD8	10-001101	IgG2 _a	Epitope has not been mapped	Unlabeled	100µg	\$210
Mouse anti T-cell receptor CD8	61-001102	IgG2 _a	Epitope has not been mapped	HRP labeled	50µg	\$310
Mouse anti T-cell receptor CD8	61-001103	IgG2 _a	Epitope has not been mapped	FITC labeled	50µg	\$310
Mouse anti T-cell receptor CD8	61-001105	IgG2 _a	Epitope has not been mapped	Biotin labeled	50µg	\$310

Measles Virus Monoclonal Antibodies

Antibody	Catalog #	Isotype	Recognises:	Label	Qty	Price
Mouse anti-Measles Virus Nucleoprotein Clone: 19G1	10-1601	IgG ₁	Recognises Catalog #071001	Unlabeled	100µg	\$300
Mouse anti-Measles Virus Nucleoprotein Clone: 19G1	61-1603	IgG ₁	Recognises Catalog #071001	FITC labeled	50µg	\$395
Mouse anti-Measles Virus Nucleoprotein Clone: 19G1	61-1605	IgG ₁	Recognises Catalog #071001	Biotin labeled	50µg	\$395

HIV-1 IIIB Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Rabbit anti-gp120	20-000401	Serum	gp120 (Catalog #012001)		0.5ml	\$345
Rabbit anti-gp120	20-000402	IgG Purified	gp120 (Catalog #012001)	Unlabeled	100µg	\$286
Rabbit anti-gp120	60-00402-3	IgG Purified	gp120 (Catalog #012001)	FITC labeled	50µg	\$392
Rabbit anti-gp120	60-00402-5	IgG Purified	gp120 (Catalog #012001)	Biotin labeled	50µg	\$392
Rabbit anti-gp160	20-000411	Serum	gp160 (Catalog #024101)		0.5ml	\$445
Rabbit anti-gp160	60-0412-3	IgG Purified	gp160 (Catalog #024101)	FITC labeled	50µg	\$445
Goat anti-p17	20-000321	Serum	p17 (Catalog #049001)		0.5ml	\$313
Goat anti-p17	20-000322	IgG Purified	p17 (Catalog #049001)	Unlabeled	100µg	\$286
Goat anti-p17	60-00322-2	IgG Purified	p17 (Catalog #049001)	HRP labeled	50µg	\$366
Goat anti-p17	60-00322-3	IgG Purified	p17 (Catalog #049001)	FITC labeled	50µg	\$366
Goat anti-p17	60-00322-5	IgG Purified	p17 (Catalog #049001)	Biotin labeled	50µg	\$366
Rabbit anti-p24	20-000201	Serum	p24 (Catalog #014001)		0.5ml	\$345
Rabbit anti-p24	20-000202	IgG Purified	p24 (Catalog #014001)	Unlabeled	100µg	\$286
Rabbit anti-p24	60-00202-3	IgG Purified	p24 (Catalog #014001)	FITC labeled	50µg	\$366
Rabbit anti-p24	60-00202-5	IgG Purified	p24 (Catalog #014001)	Biotin labeled	50µg	\$366
Goat anti-p24	20-000421	Serum	p24 (Catalog #014001)		0.5ml	\$370
Goat anti-p24	20-000422	IgG Purified	p24 (Catalog #014001)	Unlabeled	100µg	\$286
Goat anti-p24	60-00422-3	IgG Purified	p24 (Catalog #014001)	FITC labeled	50µg	\$366
Goat anti-p24	60-00422-5	IgG Purified	p24 (Catalog #014001)	Biotin labeled	50µg	\$366
Sheep anti-p24	20-000701	Serum	p24 (Catalog #014001)		0.5ml	\$345
Sheep anti-p24	20-000702	IgG Purified	p24 (Catalog #014001)	Unlabeled	100µg	\$266
Sheep anti-p24	60-00702-2	IgG Purified	p24 (Catalog #014001)	HRP labeled	50µg	\$345
Sheep anti-p24	60-00702-3	IgG Purified	p24 (Catalog #014001)	FITC labeled	50µg	\$366
Sheep anti-p24	60-00702-5	IgG Purified	p24 (Catalog #014001)	Biotin labeled	50µg	\$366
Sheep anti-p24	20-000703	IgG Purified	p24 (Catalog #014001)	Peptide Mono-Specific	100µg	\$366
Rabbit anti-rev	20-000501	Serum	Rev (Catalog #016001)		0.5ml	\$455
Rabbit anti-rev	20-000502	IgG Purified	Rev (Catalog #016001)	Unlabeled	100µg	\$402
Rabbit anti-rev	60-00502-3	IgG Purified	Rev (Catalog #016001)	FITC labeled	50µg	\$366
Rabbit anti-rev	60-00502-5	IgG Purified	Rev (Catalog #016001)	Biotin labeled	50µg	\$366
Rabbit anti-tat	20-000601	Serum	tat (Catalog #017001)		0.5ml	\$445
Rabbit anti-tat	20-000602	IgG Purified	tat (Catalog #017001)	Unlabeled	100µg	\$403
Rabbit anti-tat	60-00602-3	IgG Purified	tat (Catalog #017001)	FITC labeled	50µg	\$472
Rabbit anti-tat	60-00602-5	IgG Purified	tat (Catalog #017001)	Biotin labeled	50µg	\$472
Rabbit anti-RT	20-000511	Serum	RT (Catalog #026001)		0.5ml	\$313
Rabbit anti-RT	20-000512	IgG Purified	RT (Catalog #026001)	Unlabeled	100µg	\$260
Rabbit anti-RT	60-00512-3	IgG Purified	RT (Catalog #026001)	FITC labeled	50µg	\$366

HIV-1 IIIB Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Rabbit anti-RT	60-00512-5	IgG Purified	RT (Catalog #026001)	Biotin labeled	50µg	\$366
Rabbit anti-Protease	20-000801	Serum	Protease (Catalog #031001)		0.5ml	\$295
Rabbit anti-Protease	20-000802	IgG Purified	Protease (Catalog #031001)	Unlabeled	100µg	\$245
Rabbit anti-Protease	60-00802-3	IgG Purified	Protease (Catalog #031001)	FITC labeled	50µg	\$345
Rabbit anti-Protease	60-00802-5	IgG Purified	Protease (Catalog #031001)	Biotin labeled	50µg	\$345

SIV MAC251 Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Rabbit anti-gp120	20-000111	Serum	gp120 (Catalog #033001)		0.5ml	\$345
Rabbit anti-gp120	20-000112	IgG Purified	gp120 (Catalog #033001)	Unlabeled	100µg	\$270
Rabbit anti-gp120	60-00112-3	IgG Purified	gp120 (Catalog #033001)	FITC labeled	50µg	\$270
Rabbit anti-gp120	60-00112-5	IgG Purified	gp120 (Catalog #033001)	Biotin labeled	50µg	\$350
Rabbit anti-p28	20-000211	Serum	p28 (Catalog #022001)		0.5ml	\$345
Rabbit anti-p28	20-000212	IgG Purified	p28 (Catalog #022001)	Unlabeled	100µg	\$266
Rabbit anti-p28	60-00212-3	IgG Purified	p28 (Catalog #022001)	FITC labeled	50µg	\$370
Rabbit anti-p28	60-00212-5	IgG Purified	p28 (Catalog #022001)	Biotin labeled	50µg	\$370

HIV-2 ROD Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Rabbit anti-gp105	20-000301	Serum	gp105 (Catalog #015001)		0.5ml	\$313
Rabbit anti-gp105	20-000302	IgG Purified	gp105 (Catalog #015001)	Unlabeled	100µg	\$260
Rabbit anti-gp105	60-00302-5	IgG Purified	gp105 (Catalog #015001)	Biotin labeled	50µg	\$366
Rabbit anti-p26	20-000311	Serum	p26 (Catalog #018001)		0.5ml	\$345
Rabbit anti-p26	20-000312	IgG Purified	p26 (Catalog #018001)	Unlabeled	100µg	\$286
Rabbit anti-p26	60-00312-3	IgG Purified	p26 (Catalog #018001)	FITC labeled	50µg	\$366
Rabbit anti-p26	60-00312-5	IgG Purified	p26 (Catalog #018001)	Biotin labeled	50µg	\$366
Rabbit anti-rev	20-000611	Serum	rev (Catalog #019001)		0.5ml	\$445
Rabbit anti-rev	20-000612	IgG Purified	rev (Catalog #019001)	Unlabeled	100µg	\$402
Rabbit anti-rev	60-00612-3	IgG Purified	rev (Catalog #019001)	FITC labeled	50µg	\$345
Rabbit anti-rev	60-00612-5	IgG Purified	rev (Catalog #019001)	Biotin labeled	50µg	\$345

T-Cell Receptor Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Rabbit anti-CD4	20-000101	Serum	CD4 (Catalog #013001)		0.5ml	\$345
Rabbit anti-CD4	20-000102	IgG Purified	CD4 (Catalog #013001)	Unlabeled	100µg	\$286
Rabbit anti-CD4	60-00102-3	IgG Purified	CD4 (Catalog #013001)	FITC labeled	50µg	\$366

Regulatory Factor Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Rabbit anti-NF- κ B	20-000911	Serum	Recombinant NF- κ B		0.5ml	\$345
Rabbit anti-NF- κ B	60-00912-3	IgG Purified	Recombinant NF- κ B	FITC labeled	50 μ g	\$366

Measles Virus Polyclonal Antibodies

Antibody	Catalog #	Format	Immunogen:	Label	Qty	Price
Goat anti-Measles Virus Nucleoprotein	20-000221	Serum	Catalog #071001		0.5ml	\$313
Goat anti-Measles Virus Nucleoprotein	20-000222	IgG Purified	Catalog #071001	Unlabeled	100 μ g	\$286
Goat anti-Measles Virus Nucleoprotein	60-00222-3	IgG Purified	Catalog #071001	FITC labeled	50 μ g	\$360