

*Antibodies and Antigens
are in Our Blood.*



Fitzgerald Industries International
34 Junction Square Drive
Concord MA 01742-3049 USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
Email: antibodies@fitzgerald-fii.com

Polyclonal Antibody Data Sheet

Product: 5-HT5A Receptor (HTR5A), Human

Catalog #: 20R-HR012 **Host :** Rabbit

Immunogen: Synthetic peptide - KLH conjugated

Lot #: TBD **Affinity Constant:** N/A

Titer/Concentration: 50 µg per vial

Cross Reactivity: Specific for 3rd cytoplasmic domain of human.

Form: Purified by Immunoaffinity Chromatography in PBS containing 0.1% NaN₃.

Storage: Short term 4°C, Long term -70°C

Additional Remarks: Recommended use for Immunohistochemistry: (10mg/ml)(Optimal dilution to be determined by the researcher).

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5