



Antibodies and Antigens
are in Our Blood.

Fitzgerald Industries International, Inc.
34 Junction Square Drive
Concord MA 01742-3049 USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
E-mail: antibodies@fitzgerald-fii.com

Polyclonal Antibody Data Sheet

Product: Chicken anti **6 HIS**

Catalog #: RDI-PHIS1abCK

New Catalog #: 20R-HC001

Host: Chicken

Immunogen: Chicken were immunized with 6-HIS (HHHHHH) conjugated to KLH

Affinity Constant: Not tested

Titer/Concentration: 1 mg/ml
ELISA ; Coating 1:100-1:500 Primary 1:1,000-1:30,000
ImmunoCytochemistry : 1:100- 1:400
Western Blot 1:1,000 - 1:30,000

Cross Reactivity:

Form: Liquid in PBS,pH 7.2 with 0.1% Sodium Azide

Storage: 4°C

Additional Remarks: Suitable for Western Blot, ELISA and Histochemistry

FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.