



Antibodies and Antigens
are in Our Blood.

Fitzgerald Industries International, Inc.
34 Junction Square Drive
Concord MA 01742-3049 USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
E-mail: antibodies@fitzgerald-fii.com

Polyclonal Antibody Data Sheet

Product: Goat anti **8-Hydroxyguanosine**

Catalog #: RDI-8OHGUANabG

New Catalog #: 20R-HG001 **Host:** Goat

Immunogen: 8-Hydroxyguanosine-conjugate

Affinity Constant:

Titer/Concentration: 0.25 ml per Vial whole serum

Cross Reactivity: Amino Acid sequence 42-59 of human receptor for advanced glycated end products

Form: Whole serum with no preservative

Storage: 4°C Short term. Long term store at -20°C

Additional Remarks: ELISA 1:250,000
IHC using formalin-fixed, paraffin-embedded sections ,1:200

FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.