

*Antibodies and Antigens
are in Our Blood.*

Fitzgerald Industries International
Suite 1A North, 30 Sudbury Road
Acton, MA 01720, USA
Tel: 978.371.6446
Fax: 978.371.2266



www.fitzgerald-fii.com
Email: antibodies@fitzgerald-fii.com

Normal Animal Serum

Product: Red Blood Cells, Chicken Fresh 10% Conc. in Alsever's

Catalog #: 88R-C001 **Grade :**

Lot #: TBD **Purity:**

Source: Fresh Chicken Red Blood Cells

Concentration: 25 ml per vial

Contaminants: The whole blood is aseptically collected and prepared.

Form: Fresh Chicken red blood cells stabilized 10% Concentration in Alsever's.

Biohazard: Not for human use, use normal laboratory precautions when handling reagent. For research purposes only.

Storage: Store at 4°C.

Additional Remarks: The whole blood is aseptically collected and prepared. Washed cell suspensions are useful for the titration of complement, adsorption procedures, testing for agglutinins, and for the preparation of stroma as particulate reagents.

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user names assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5