

*Antibodies and Antigens  
are in Our Blood.*

Fitzgerald Industries International  
Suite 1A North, 30 Sudbury Road  
Acton, MA 01720, USA  
Tel: 978.371.6446  
Fax: 978.371.2266



www.fitzgerald-fii.com  
Email: antibodies@fitzgerald-fii.com

## **Normal Animal Serum**

**Product:** Red Blood Cells, Chicken, Glut Stabilized

**Catalog #:** 88R-C002 **Grade :**

**Lot #:** TBD **Purity:**

**Source:** Chicken Red Blood Cells

**Concentration:** 20 ml per vial

### **Contaminants:**

**Form:** Chicken Red Blood Cells suspended in saline with 0.1 % Sodium Azide . Cells are produced from freshly washed Chicken RBCs by brief exposure to glutataldehyde-buffer solution and exhaustive washing in saline solution.

**Biohazard:** Use standard laboratory precautions when handling this material.

**Storage:** Store at 4°C .

**Additional Remarks:** Treatment largely preserves the native antigenicity of the erythrocytes while rendering the cells generally resistant to lysis by osmotic shock, freeze-thawing or immune hemolysis. These cells are used directly in hemagglutination procedures, or may be coupled with various proteins for hemagglutination testing.

#### **FOR INVITRO USE ONLY**

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user names assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5