

*Antibodies and Antigens  
are in Our Blood.*



Fitzgerald Industries International  
Suite 1A North, 30 Sudbury Road  
Acton, MA 01720, USA  
Tel: 978.371.6446  
Fax: 978.371.2266

www.fitzgerald-fii.com  
Email: antibodies@fitzgerald-fii.com

## Normal Animal Serum

**Product:** Normal Equine Serum

**Catalog #:** 88-NE21 **Grade :** N/A

**Lot #:** S10011903 **Purity:** N/A

**Source:** Normal Equine Blood

**Concentration:** Total Protein: 6.6 g/dl Hemoglobin: 6.8 mg/dl

**Contaminants:** Endotoxin: <0.125 EU/ml (Limulus Amebocyte Lysate Assay)  
Mycoplasma: None Detected (Modified Barile Method)  
Ba

**Form:** 0.2µm filtered serum, pH 7.5

**Biohazard:** Use standard Laboratory procedures.

**Storage:** 4°C for shipping purposes, for routine storage aliquot and freeze at -20°C or below

**Additional Remarks:** Produced from quarantined donor horses that are tested in accordance with USDA Veterinary Services memorandum #800.71 for Equine Infectious Anemia (Coggins test).

for in vitro laboratory use only

**FOR INVITRO USE ONLY**

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user names assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5