



## Antigen Data Sheet

<b><u>Product:</u></b>	Factor I Inactivator		
<b><u>Catalog #:</u></b>	30C-CP6028	<b><u>Grade :</u></b>	Highly pure
<b><u>Lot #:</u></b>	TBD	<b><u>Purity:</u></b>	>95% by SDS PAGE
<b><u>Source:</u></b>	Prepared from human serum/plasma		
<b><u>Concentration:</u></b>	batch dependent		
<b><u>Contaminants:</u></b>	Only trace amounts of IgG, IgA, IgM, Factor H, Factor B, Factor P, C3, C4, or albumin are detected.		
<b><u>Form:</u></b>	Supplied as a liquid in PBS, pH 7.4. NO PRESERVATIVE ADDED.		
<b><u>Biohazard:</u></b>	This product has been tested at serum donor level by FDA approved methods and found NEGATIVE for HIV I/II, and HCV antibodies and Hepatitis B surface antigen. It is recommended to observe full safety precautions when handling this material, since no testing method can completely assure the absence of all infectious agents.		
<b><u>Storage:</u></b>	-70°C		
<b><u>Additional Remarks:</u></b>	280nm $\epsilon$ = 1.43 1cm, 0.1%		

**FOR INVITRO USE ONLY**

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5