



Antigen Data Sheet

Product: *Factor IX, Human, Antigen*

Catalog #: RDI-FAC9-AG

New Catalog #: 30R-AH037

Grade: Highly purified

Purity: >95 % by SDS PAGE

Source: Human Plasma

Concentration: 1.56 mg/ml

Contaminants:

Form: Supplied in 20 mM Tris-HCL, 0.1 M NaCl, 1 mM Benzamidine pH 7.4

Biohazard: Starting material has tested negative for HBsAg, HCV and HIV 1/2 by FDA approved test methods. Donors are screened for CJD (Creutzfeldt-Jakob Disease)

Storage: -70°C ,

Additional Remarks: Activity 185.00 PEU/ mg

FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.