



## Antigen Data Sheet

**Product:** Factor V, Bovine

**Catalog #:** 30C-CP3105U      **Grade :** Highly pure

**Lot #:** TBD      **Purity:** >98% by SDS PAGE

**Source:** Prepared from bovine serum/plasma

**Concentration:** batch dependent

### **Contaminants:**

**Form:** Supplied as a liquid in 50% (vol/vol) glycerol/deionized water containing 2 mM Calcium Chloride. NO PRESERVATIVE ADDED.

**Biohazard:** This product is produced from non-human material. As with any material produced from animal tissue, it is recommended to observe full safety precautions when handling this material.

**Storage:** -20°C

**Additional Remarks:** Molecular Weight: 330,000 daltons

280nm  
 $\epsilon = 0.96$   
1cm, 0.1%

### **FOR INVITRO USE ONLY**

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5