



## Antigen Data Sheet

**Product:** *Factor VIIA, Human*

**Catalog #:** RDI-FAC7A-AG

**New Catalog #:** 30R-AF029

**Grade:** Highly purified

**Purity:** >95 % by SDS PAGE

**Source:** Human Plasma

**Concentration:** 1.48 mg/ml

**Contaminants:**

**Form:** Supplied in 20 mM Tris-HCL, 0.1 M NaCl, pH 7.4

**Biohazard:** Starting material has tested negative for HBsAg, HCV, HIV 1/2, Anti-HTLV I&II, and Syphilis by FDA approved test methods. Donors are screened for CJD (Creutzfeldt-Jakob Disease)

**Storage:** -70°C,

**Additional Remarks:** Activity 32189 PEU/ mg

**FOR INVITRO USE ONLY**

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.