



Antigen Data Sheet

Product: *Bovine Factor X*

Catalog #: CP2100U

New Catalog #: 30C-CP2100U **Grade:** Ultra Pure

Lot #: A08121750 **Purity:** > 98% by SDS - PAGE

Source: Prepared from bovine serum / plasma

Concentration: 8.8 mg/ml

Contaminants:

Form: Supplied as a liquid in 50% (vol / vol) glycerol / deionized water.

Biohazard: Not for human use. Use normal laboratory precautions when handling reagent. For in vitro research purposes only.

Storage: Store at -20°C.

Additional Remarks:

FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.