



Antibodies and Antigens
are in Our Blood.

Fitzgerald Industries International, Inc.
34 Junction Square Drive
Concord MA 01742-3049 USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
E-mail: antibodies@fitzgerald-fii.com

Normal Animal Serum

Product: *Fetal Bovine Serum*

Catalog #: 88-NB05

New Catalog #: 88-NB05

Grade: N/A

Purity: N/A

Source: Fetal Bovine Serum

Concentration: N/A

Contaminants: Endotoxin: 0.15 ng/ml
Negative for Mycoplasma, Bacteria and Fungi
Viral Agents tested: BVD, IBR and PI-3 not detected

Form: Sterile filtered. pH 7.3 BIOCHEMICAL ANALYSIS: AST 23.0 IU/L; ALT 8.0 IU/L;
Alkaline Phosphatase 211.0 IU/L; Bilirubin 0.4 mg/dL; Protein 4.0 g/dL;
Albumin 3.0 g/dL; Na 122.0 mEq/L; K 17.0 mEq/dL; P 2.42 mg/dL;
Glucose 29.0 mg/dL; BUN 10.7 mg/dL; Creatinine 1.76 mg/dL; IgG 63 ug/ml;
Cholesterol 14mg/dL; Osmolality 305 mOsm/Kg; Hemoglobin 5.9 mg/dL

Biohazard: Use standard Laboratory procedures.

Storage: -20°C or below

Additional Remarks:

FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.