



Polyclonal Antibody Data Sheet

Product: Acid Phosphatase

Catalog #: 20C-CR6098RP **Host :** Rabbit

Immunogen: Acid Phosphatase [Potato]

Lot #: TBD **Affinity Constant:** N/A

Titer/Concentration: 1ml per vial

Cross Reactivity: Assay by immunoelectrophoresis resulted in a single precipitin arc against anti-Rabbit Serum as well as purified and partially purified Acid Phosphatase [Potato]. Cross reactivity against Acid Phosphatase from other tissues and species may occur but have not been specifically determined.

Form: Supplied in 0.02 M Potassium Phosphate, 0.15 M Sodium Chloride, pH 7.2, and 0.01% Sodium Azide as a preservative.

Storage: 4°C short term, -20°C long term

Additional Remarks: ELISA 1:4,000 - 1:40,000
WESTERN BLOT 1:2,000 - 1:10,000
IMMUNOHISTOCHEMISTRY 1:1,000 - 1:5,000

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.