

*Antibodies and Antigens
are in Our Blood.*

Fitzgerald Industries International
Suite 1A North, 30 Sudbury Road
Acton, MA 01720, USA
Tel: 978.371.6446
Fax: 978.371.2266



www.fitzgerald-fii.com
Email: antibodies@fitzgerald-fii.com

Polyclonal Antibody Data Sheet

Product: Actin, Chicken/Rat

Catalog #: 20C-CR2001R **Host :** rabbit

Immunogen: Chicken back muscle actin bound to KLH.

Lot #: TBD **Affinity Constant:** N/A

Titer/Concentration: batch dependent

Cross Reactivity: The antiserum has been evaluated for its ability to give a specific and positive immunostaining of actin-containing cells in cultures of rat fibroblasts using both the immunoperoxidase and immunofluorescence procedures.

Form: Supplied as a liquid containing 0.1% Sodium Azide

Storage: -20°C

Additional Remarks:

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5