



Antibodies and Antigens  
are in Our Blood.

Fitzgerald Industries International, Inc.  
34 Junction Square Drive  
Concord MA 01742-3049 USA  
Tel: 978.371.6446  
Fax: 978.371.2266

www.fitzgerald-fii.com  
E-mail: antibodies@fitzgerald-fii.com

## Normal Animal Serum

<b><u>Product:</u></b>	<b><i>Normal Mouse Serum</i></b>		
<b><u>Catalog #:</u></b>	88-NM35	<b><u>Grade:</u></b>	Neat serum
<b><u>New Catalog #:</u></b>	88-NM35	<b><u>Purity:</u></b>	N/A
<b><u>Source:</u></b>	Normal mice		
<b><u>Concentration:</u></b>	Total protein 5.8 g/dl pH 8.1		
<b><u>Contaminants:</u></b>	Hemoglobin 17.29mg/dl		
<b><u>Form:</u></b>	0.2 $\mu$ m sterile filtered. No preservatives added.		
<b><u>Biohazard:</u></b>	Use normal laboratory safety techniques when handling this product.		
<b><u>Storage:</u></b>	-20°C for best results		
<b><u>Additional Remarks:</u></b>			

### FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.