

*Antibodies and Antigens
are in Our Blood.*



Fitzgerald Industries International
Suite 1A North, 30 Sudbury Road
Acton, MA 01720, USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
Email: antibodies@fitzgerald-fii.com

Normal Animal Serum

Product: Bovine Red Blood Cells, Glut Stabilized

Catalog #: 88R-B002 **Grade :**

Lot #: TBD **Purity:**

Source: Bovine Red Blood Cells

Concentration: 20 ml/vial

Contaminants:

Form: Bovine Red Blood Cells suspended in Saline with 0.1 % Sodium Azide . Cells are produced from freshly washed Bovine RBC by brief exposure to glutaraldehyde-buffer solution and exhaustive washing in saline solution.

Biohazard: Not for human use, use normal laboratory precautions when handling reagent. For research purposes only.

Storage: 4°C

Additional Remarks: Brief exposure to a glutaraldehyde-buffer solution and exhaustive washing in saline largely preserves the native antigenicity of the erythrocytes while rendering the cells generally resistant to lysis by osmotic shock, freeze-thawing or immune hemolysis.

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user names assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5