



Antibodies and Antigens
are in Our Blood.

Fitzgerald Industries International, Inc.
34 Junction Square Drive
Concord MA 01742-3049 USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
E-mail: antibodies@fitzgerald-fii.com

Normal Animal Serum

<u>Product:</u>	<i>Serum, Rabbit, IgG, Rat cell adsorbed</i>		
<u>Catalog #:</u>	RDI-SDG403-R		
<u>New Catalog #:</u>	89R-RR001	<u>Grade:</u>	IgG Fraction
		<u>Purity:</u>	
<u>Source:</u>	Serum from Normal Non- immunized Rabbits processes from whole blood.		
<u>Concentration:</u>	5 mg/ml		
<u>Contaminants:</u>			
<u>Form:</u>	Normal rabbit serum adsorbed with Rat red blood cells, thymocytes and spleen cells. IgG Fractionated then 0.2 μ filtered		
<u>Biohazard:</u>	Use normal laboratory practices and procedures.		
<u>Storage:</u>	Divide into useful aliquots and Store at -20°C		
<u>Additional Remarks:</u>	Remove under sterile conditions or add preservative.		

FOR INVITRO USE ONLY

CAUTION, not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness of purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.