

*Antibodies and Antigens
are in Our Blood.*



Fitzgerald Industries International
Suite 1A North, 30 Sudbury Road
Acton, MA 01720, USA
Tel: 978.371.6446
Fax: 978.371.2266

www.fitzgerald-fii.com
Email: antibodies@fitzgerald-fii.com

Normal Animal Serum

Product: Rat Fresh Macrophages

Catalog #: 88R-R006

Grade :

Lot #: TBD

Purity: Viability by trypan blue exclusion > 99%

Source: Sprague-Dawley derived rats

Concentration: 2 x 10⁷ washed viable cells

Contaminants:

Form: Fresh rat macrophages in RPMI-1640 with 100mcg neomycin per ml

Biohazard: Not for human use, use normal laboratory precautions when handling reagent. For research purposes only.

Storage: Store at 4°C.

Additional Remarks: These cells are useful to prepare lymphokines and as starting material to isolate lymphocyte subsets

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user names assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5