



## Highly Purified Polyclonal Antibody

### Data Sheet

**Product:** 4-1BB Receptor

**Catalog #:** 70R-BG004      **Host :** goat

**Immunogen:** highly pure (>98%) recombinant human 4-1BB Receptor

**Lot #:** TBD      **Affinity Constant:** N/A

**Titer/Concentration:** 50 µg/vial; lyophilized

**Cross Reactivity:**

**Form:** lyophilized from PBS, pH 7.2

**Storage:** 4°C short term, -20°C long term

**Additional Remarks:** purified by affinity chromatography employing immobilized h4-1BB Receptor matrix

**FOR INVITRO USE ONLY**

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.