



Highly Purified Polyclonal Antibody

Data Sheet

Product: 4-1BBL

Catalog #: 70R-BR007 **Host :** rabbit

Immunogen: highly pure (>98%) recombinant human 4-1BBL

Lot #: TBD **Affinity Constant:** N/A

Titer/Concentration: 50 µg/vial; lyophilized

Cross Reactivity: specific for anti-human 4-1BBL

Form: lyophilized from PBS, pH 7.2

Storage: 4°C short term; -20°C long term

Additional Remarks: purified by affinity chromatography employing immobilized h4-1BBL matrix

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.