



Polyclonal Antibody Data Sheet

<u>Product:</u>	5-HT1F Receptor (HTR1F), Human		
<u>Catalog #:</u>	20R-HR015	<u>Host :</u>	Rabbit
<u>Immunogen:</u>	Synthetic peptide - KLH conjugated		
<u>Lot #:</u>	TBD	<u>Affinity Constant:</u>	N/A
<u>Titer/Concentration:</u>	50 µg per vial		
<u>Cross Reactivity:</u>	Specific for 3rd cytoplasmic domain of human.		
<u>Form:</u>	Purified by Immunoaffinity Chromatography in PBS containing 0.1% NaN ₃ .		
<u>Storage:</u>	Short term 4°C, Long term -70°C		
<u>Additional Remarks:</u>	Recommended use for Immunohistochemistry: (6.5mg/ml)(Optimal dilution to be determined by the researcher).		

FOR INVITRO USE ONLY

CAUTION: Not for use in humans or clinical diagnosis: This product is intended for research or manufacturing use only. It is pharmaceutically unrefined, and verification of its suitability for use in humans or as clinical diagnostic reagents and the compliance with all Federal and State laws regulating such applications are the responsibility of the purchaser. There is no expressed or implied warranty. No liability is assumed for fitness or purpose and merchantability, or direct or consequential damages. The user assumes all responsibility for care, custody and control of the material, including its disposal, in accordance with all regulations.

FM 4.15-5