

NEW



ZIKA VIRUS PRODUCTS

Monoclonal Zika NS1 and Envelope Antibodies

Polyclonal Zika NS1 Antibodies

Recombinant Zika NS1 and Envelope Proteins*

*Expressed in Insect Cells & E.coli

- ◇ Tested in ELISA and Lateral Flow
- ◇ Free Samples & Bulk Pricing Available

Check out www.fitzgerald-fii.com and search for: **Zika**



Fitzgerald

www.fitzgerald-fii.com
T 978.371.6446 F 978.371.2266
antibodies@fitzgerald-fii.com

MONOCLONAL ANTIBODIES

Catalog No.	Product Name	Host	Clone	Applications	URL	Price (1 mg)	Price (500 ug)
10-2707	Zika virus NS1 antibody	Mouse	M86640	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2708	Zika virus NS1 antibody	Mouse	M86641	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2709	Zika virus NS1 antibody	Mouse	M86642	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2710	Zika virus NS1 antibody	Mouse	M86643	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2711	Zika virus NS1 antibody	Mouse	M86644	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2712	Zika virus envelope antibody	Mouse	M86645	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2713	Zika virus envelope antibody	Mouse	M86646	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2714	Zika virus envelope antibody	Mouse	M86647	ELISA, Lateral Flow	Click here	\$495.00	\$295.00
10-2715	Zika virus envelope antibody	Mouse	M86648	ELISA, Lateral Flow	Click here	\$495.00	\$295.00

POLYCLONAL ANTIBODIES

Catalog No.	Product Name	Host	Applications	URL	Price (500 ug)
70-1086	Zika virus NS1 antibody	Goat	ELISA	Click here	\$225.00
70-1087	Zika virus NS1 antibody	Goat	ELISA	Click here	\$425.00

RECOMBINANT PROTEINS

Catalog No.	Product Name	Expression System	Purity	Applications	URL	Price (1 mg)	Price (100 ug)
30-1932	Zika virus envelope protein	Insect cells	> 85% pure	ELISA, Lateral Flow	Click here	\$2,995.00	\$325.00
30-1935	Zika virus envelope protein	E.coli	> 95% pure	ELISA, WB, Lateral	Click here	\$1,595.00	\$195.00
30-1931	Zika virus NS1 protein	Insect cells	> 90% pure	ELISA, Lateral Flow	Click here	\$2,995.00	\$325.00
30-1936	Zika virus NS1 protein	E.coli	> 95% pure	ELISA, WB, Lateral	Click here	\$1,595.00	\$195.00

MATCHED PAIR INFORMATION

- Antibodies 10-2712, 10-2713, 10-2714 & 10-2715 can be used in assays to detect Protein 30-1932
- Antibodies 10-2707, 10-2708, 10-2709 & 10-2710 can be used in assays to detect Protein 30-1931